EN Instructions for Use

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Description

Compoglass® F is a light-curing, radiopaque compomer-based restorative material. Compoglass F cures with light in the wavelength range of 400–500 nm (blue light). Compoglass F combines the advantages of glass ionomer cements with those of light-curing composite materials.

Shades

Compoglass F is available in several shades.

Composition

The monomer matrix is composed of dimethacrylates (19.30 wt.%). The inorganic fillers include ytterbium trifluoride, Ba-Al-fluorosilicate glass and spheroid mixed oxide (80.5 wt.%). Additionally, initiators, stabilizers and pigments (0.20%) are included. The total content of inorganic fillers is 57.5 Vol.% or 80.5 wt.% respectively. The particle size is between 0.2 and 3.0 µm.

Indication

- Restoration of deciduous teeth
- Class V restorations (cervical caries, root erosion, wedge-shaped defects)
- Class III restorations
- small to medium sized restorations of classes I and II

Contraindication

The placement of Compoglass F restorations is contraindicated

- if a patient is known to be allergic to any of the ingredients of Compoglass F;
- if the stipulated working technique cannot be applied;
- for direct and indirect pulp capping;
- for permanent, stress-bearing occlusal restorations in teeth of the secondary dentition.

Side effects

In rare cases, ingredients of Compoglass F may cause a sensitizing reaction. The product should not be used in such cases. To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (apply a calcium hydroxide-based preparation to areas in close proximity of the pulp) and covered with an adequate lining.

Interactions

Substances containing eugenol/clove oil may inhibit the polymerization of compomers. These materials must not be used in conjunction with Compoglass F. Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

CE 0123

Rx ONLY

Date information prepared:

2014-11-10/Rev. 3

573127/WE3

Manufacturer
Ivoclar Vivadent AG
FL-9494 Schaan/Liechtenstein
www.ivoclarvivadent.com

ivoclar
vivadent
clinical

Usage

1. Shade determination

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist.

2. Isolation

Appropriate isolation is required, preferably with a rubber dam (e.g. OptraDam® Plus).

3. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e., by preserving as much of the tooth structure as possible. Do not prepare sharp internal edges or angles or additional undercuts in caries-free areas. The dimensions of the cavity are generally dictated by the extent of the caries or the size of the previous restorations. Bevel enamel edges of anterior teeth. In the posterior region, only the sharp enamel edges should be rounded. Do not prepare caries-free cervical defects. Clean them with pumice or a suitable cleansing paste (e.g. Proxylt® RDA 36, medium) and with the help of rubber polishers or rotary brushes. Ablate sclerotised dentin around the tooth neck in a thin layer with a rotary instrument (diamond finisher, coarse disc). Subsequently, rinse the cavity with water to remove all the residue and dry it with water- and oil-free air.

4. Pulp protection / Base

Do not apply a base if an enamel-dentin bonding agent is used. In very deep cavities, areas close to the pulp should be selectively coated with a calcium hydroxide liner (e.g. ApexCal®) and subsequently covered with pressure-resistant cement (glass ionomer cement, e.g. Vivaglass® Liner). Do not cover the remaining cavity walls so that they can be used to generate a bond with an enamel-dentin adhesive.

5. Apply matrix / interdental wedge

Use a wrap around matrix for cavities affecting the proximal area or a sectional matrix and wedge it.

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the instructions for use of the product in use. Ivoclar Vivadent recommends using the adhesive Adhese® Universal, for example.

7. Application of Compoglass F

- In order to achieve optimum results, Compoglass F should be applied in layers of max. 3 mm (light shades) or 2 mm (dark shade = 340/A4) and adapted with the suitable instruments (e.g. OptraSculpt®). Use a special contact point instrument (e.g. OptraContact) for large cavities where contacts are difficult to create.
- For complete polymerization, sufficient exposure time is required. For the recommendations regarding exposure time and light intensity see Table 1.
- Hold the light emission window of the curing light as close to the restorative as possible.
- When using a metal matrix, additionally polymerize from the buccal or the lingual/palatal aspect after removing the matrix.

Unit	Bluephase C8	Bluephase Style	Bluephase 20i
Program			
Turbo Program	–	–	2x 5 s
High Power Program	2x 20 s	20 s	20 s
Soft Start Program	2x 20 s	–	20 s
Light intensity		Exposure time	
≥ 600 mW/cm²		40 s	
≥ 1,000 mW/cm²		20 s	

Table 1

8. Finishing / Checking the occlusion / Polishing

After polymerization, remove excess material with suitable finishers (e.g. Astropol® F) or fine diamonds. Remove proximal excess material with flexible discs, finishing strips or diamond-coated files. Check the occlusion and articulation and grind in the necessary adjustments to prevent premature contacts or undesired occlusal paths on the surface of the restoration. Use silicone polishers (e.g. OptraPol®) as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Additional notes

1. Use Compoglass F at room temperature. Cold material may be difficult to dispense.
2.  For single use only. If Compoglass F is directly applied from the Cavifill in the mouth of the patient, the Cavifill should only be used for one patient for hygienic reasons (prevention of cross-contamination between patients).
3. Do not disinfect Cavifills with oxidizing disinfectants.

Warning

Prevent any contact of uncured Compoglass F with skin / mucous membrane and eyes. Unpolymerized Compoglass F may cause slight irritation and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage and hygiene

- Do not use Compoglass F after the indicated date of expiration.
- Expiry date: see note on packaging.
- Storage temperature 2–28 °C (36–82 °F).
- Close Cavifills immediately after usage. Exposure to light causes premature polymerization.

Keep out of the reach of children!

For dental use only!

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.