GC Fuji IX GP

RADIOPAQUE POSTERIOR GLASS IONOMER RESTORATIVE CEMENT

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Class I and II restorations in deciduous teeth.
- Non-load bearing Class I and Class II restorations in permanent teeth.
- Intermediate restorative and base material for heavy stress situation in Class I and Class II cavities using sandwich laminate technique.
- Class V and root surface restorations.
- Core build-up.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced,discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

Powder / Liquid Ratio (g/g)	3.6 / 1.0
Mixing Time (sec.)	25-30"
Working Time (min., sec.)	2'00"
Net Setting Time (min., sec.)	2'20"

Test conditions : Temperature (23 +/-1°C)
Relative humidity (50 +/-5%)
ISO 9917-1 : 2003 (E) (Glass polyalkenoate cement)

1. POWDER AND LIQUID DISPENSING

- The standard powder to liquid ratio is 3.6g/1.0g. (1 level scoop of powder to 1 drop of liquid).
- For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
- Hold the liquid bottle vertically and squeeze gently.
- Close bottles tightly immediately after use.

2. MIXING

Dispense powder and liquid onto the pad. Using the plastic spatula, divide the powder into 2 equal parts. Mix the first portion with all the liquid for 10 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-20 seconds (Fig. 1).

3. RESTORATIVE TECHNIQUE

- Prepare the tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping, use calcium hydroxide.
- Apply GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) to the bonding surfaces using a cotton pellet or sponge.
- Rinse thoroughly with water. Dry by blotting with a cotton pellet or gently blowing with an air syringe. DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening) (Fig. 2).
- Mix the required amount of GC Fuji IX GP. Working time is 2 minutes from the start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
- Transfer cement to the preparation using a syringe or other suitable instrument. Avoid air bubbles.
- Form the preliminary contour and cover with a matrix if required (Fig. 3).
- When set, immediately apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure) (Fig. 4).

4. FINISHING

- Final finishing under water spray using standard techniques can begin at 6 minutes after starting the mix (Fig. 5).
- Apply a final coat of GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure for 10 sec.) or G-COAT PLUS (light cure for 20 sec.) (Fig. 6).

- Instruct the patient not to apply pressure for 1 hour.

STORAGE

Store at temperature of 4-25°C (39.2-77.0°F).
(Shelf life : 3 years from date of manufacture).

SHADES

A2, A3, A3.5, B2, B3, C4

Shade numbers according to Vita® shade guide.

©Vita is a trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

- 1-1 package : 15g powder, 8g (6.4mL) liquid, powder scoop, mixing pad (No. 22), plastic spatula.
- Bottle of 15g powder with scoop.
- Bottle of 8g (6.4mL) liquid.

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with any other glass ionomer components.

Last revised : 06/2008

CE 0086

GC

MANUFACTURED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

EU: GC EUROPE N.V.
Research-Park, Interleuvenlaan 13, B-3001 Leuven, Belgium
TEL: +32 16 39 80 50

GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
TEL: +1-708-597-0900

GC ASIA DENTAL PTE. LTD.
19 Loyang Way, #06-27 Singapore 508724
TEL: +65 6546 7588

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