GC Fuji TRIAGE®

RADIOPAQUE GLASS IONOMER PROTECTION AND TEMPORARY RESTORATIVE MATERIAL

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

1. Fissure protection.
2. Root surface protection.
3. Hypersensitivity prevention and control.
4. Protection for immature enamel.
5. Temporary filling of endodontic access.
6. Intermediate restorations.

CONTRAINDICATIONS

1. Pulp capping
2. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

PINK Shade is command set, WHITE Shade is chemical set only.

Powder / Liquid Ratio (g / g)	1.8 / 1.0
Mixing Time (sec.)	20-25"
Working Time (min., sec.)	1' 40"
Net Setting Time (min., sec.)	2' 30"
Final Finishing Commencing Time	6'00"
Final Finishing Commencing Time (if light cured (PINK Command set)	4'00"

Test conditions : Temperature (23±1°C), Relative humidity (50±1.5%)
ISO 9917-1 : 2003 (Dental water-based cements) (Restorative cements)

A.FISSURE PROTECTION ROOT SURFACE PROTECTION HYPERSENSITIVITY PREVENTION AND CONTROL PROTECTION FOR IMMATURE ENAMEL

1. Preparation of the tooth surfaces (e.g. for fissure protection or root surface protection)
 - a) After cleaning the tooth surfaces (prophylaxis with pumice and water) in usual manner, rinse thoroughly with water. Avoid aggravating the operculum.
Note :
If extra retention is desired, application of GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) is recommended.
 - b) Dry by blotting with a cotton pellet or gently blowing with an air syringe (Fig. A-1). DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening).
2. Powder and Liquid Dispensing
 - a) The standard powder to liquid ratio is 1.8g / 1.0g. (1 level scoop of powder to 1 drop of liquid).
 - b) For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
 - c) Make sure the liquid nozzle is clean and dry before dispensing liquid. Turn the liquid bottle horizontal and hold in this position briefly to remove air bubbles. Then invert and hold the liquid bottle vertically and squeeze gently to dispense a bubble free drop of liquid. After dispensing, wipe any residual liquid off the nozzle.
 - d) Close bottle caps tightly immediately after use.
3. Mixing
Dispense the required amounts of powder and liquid onto the pad. Using the plastic spatula, divide the powder into 2 equal parts. Spread the liquid over the pad (Fig. A-2) and mix the first portion with all the liquid for 10 seconds. Incorporate the remaining portion and mix the whole amount thoroughly for 10-15 seconds (total time within 25 seconds). Working time is 1 minutes 40 seconds from the start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time.
4. Placement
 - a) Brush and mix the material using a suitable placement instrument or brush and apply to the tooth surface (Fig. A-3). Then use a brush to spread a thin film of GC Fuji TRIAGE® directly over the root surface or hypersensitive area or over the occlusal surface and into the pits and fissures.
Note :
If a faster set is desired, use a visible light curing device* for 20-40 seconds. Place light source as closely as possible to the cement surface. This function applies only to the Pink Shade. After light cure, it is advisable to protect the surface with a varnish.
 - b) After placement, when the material starts to lose the glossy appearance (or after curing with the light curing device), apply GC Fuji VARNISH (blow dry) or GC Fuji COAT LC (light cure) to the sealed area and the margins using a cotton pellet or sponge (Fig. A-4).
 - c) Finishing under air water spray can be performed 6 minutes from start of mix (chemically set) or 4 minutes if light cured. Use a superfine diamond bur or a silicone finishing point (Fig. A-5).
 - d) Apply GC Fuji VARNISH or GC Fuji COAT LC to the area again (Fig. A-6).

B.TEMPORARY FILLING OF ENDODONTIC ACCESS

1. Cleaning the pulp chamber
 - a) After appropriate pulp treatment, clean and gently dry the pulp chamber with an air syringe (Fig. B-1).
 - b) Fill the chamber with a cotton pellet (Fig. B-2).
2. Powder / Liquid Dispensing and Mixing
See directions in A (above), Section 2 and 3.
3. Placement
 - a) Using a suitable instrument, place the mixed material over the cotton pellet (Fig. B-3).
Note :
If a faster set is desired, use a visible light curing device* for 20-40 seconds. Place light source as closely as possible to the cement surface. This function applies only to the Pink Shade. After light cure, it is advisable to protect the surface with a varnish.
 - b) Moisture protection and Finishing
See directions in A, Section 4, b, c and d.

C. INTERMEDIATE RESTORATIONS

1. Caries removal
 - a) Remove caries and loose debris with hand instruments.
 - b) For better retention, it is recommended to gently clean the carious surface with GC CAVITY CONDITIONER for 10 seconds or GC DENTIN CONDITIONER for 20 seconds.
 - c) Rinse thoroughly with water. Dry by blotting with a cotton pellet or gently blowing with an air syringe (Fig. C-1). DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening).
2. Powder / Liquid Dispensing and Mixing
See directions in A (above), Section 2 and 3.
3. Placement
 - a) Take the mixed material using a suitable placement instrument or brush and apply to the prepared lesion or tooth surface without incorporating air bubbles (Fig. C-2). Form the contour and if possible cover with a matrix.
Note :
If a faster set is desired, use a visible light curing device* for 20-40 seconds. Place light source as closely as possible to the cement surface. This function applies only to the Pink Shade. After light cure, it is advisable to protect the surface with a varnish.
 - b) See directions in A, Section 4, b, c and d.

* NOTE :

The initial set of GC Fuji TRIAGE® PINK can be accelerated using the energy from a dental halogen light curing device.

SHADE

Pink, White

STORAGE

Store at temperature of 4-25° (39.2-77.0°F)
(Shelf life : 3 years from date of manufacture)

PACKAGE

15 g powder, 10 g (8.0 mL) liquid, Dentin Conditioner 6 g (5.7 mL), powder scoop, mixing pad No. 22, plastic spatula.

CAUTION

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. DO NOT mix powder or liquid with any other glass ionomer components.

US Patent : 6264472 6756421

UK Patent : 2353042 2357773

France Patent : 2797396 2799954

Australia Patent : 768901 775349

Last revised : 05/2010

0086

GC

MANUFACTURED BY
GC CORPORATION
715-1 Hasegawa-cho, Isahaya-shi, Tokyo 174-8585, Japan
DISTRIBUTED BY
GC CORPORATION
715-1 Hasegawa-cho, Isahaya-shi, Tokyo 174-8585, Japan

EU, GC EUROPE N.V.
Research Park, Herengracht 13, B-3001 Leuven, Belgium
Tel. : +32 (0) 3 235 3042
GC AMERICA INC.
3727 West 127th Street, Anso, IL 60803 U.S.A.
Tel. : 1-708-697-8800

GC ASIA DENTAL PTE. LTD.
10 Loyang Way, #06-27 Singapore 508724
Tel. : +65 5446 7588

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