Description

Tetric EvoCeram is a light-curing, radiopaque nano-hybrid composite for direct restorations. Tetric EvoCeram cures with light in the wavelength range of 400–500 nm.

Shades

Tetric EvoCeram is available in the following 22 shades:

- Enamel shades: A1, A2, A3, A3.5, A4, B1, B2, B3, C1, C2, C3, D2, D3, T
- Dentin shades: A2 Dentin, A3.5 Dentin, A4 Dentin, B2 Dentin
- Bleach shades: Bleach XL, Bleach L, Bleach I, Bleach M

Composition

The monomer matrix is composed of dimethacrylates (17–18 wt.%). The fillers contain barium glass, ytterbium trifluoride, mixed oxide and copolymers (82–83 wt.%). Additives, initiators, stabilizers and pigments are additional ingredients (<1.0 wt.%). The total content of inorganic fillers for the Bleach shades is 79–80 wt.% or 60–61 vol.% , for all the other shades 75–76 wt.% or 53–55 vol.%. The particle sizes of the inorganic fillers range between 40 nm and 3 µm.

Indication

- Restoration of deciduous teeth
- Restorations in the posterior region (Classes I and II)
- Anterior restorations (Classes III, IV)
- Class V restorations (cervical caries, root erosion, wedge-shaped defects)
- Veneering of discoloured anterior teeth
- Splinting of mobile teeth
- Extended fissure sealing in molars and premolars
- Repair of composite/ceramic veneers
- Build-ups for transparent, removable, Invisalign® orthodontic retainers

Contraindication

The placement of Tetric EvoCeram restorations is contraindicated:

- If a dry working field cannot be established, or the stipulated working technique cannot be applied.
- If a patient is known to be allergic to any of the ingredients of Tetric EvoCeram.

Side effects

In individual cases, components of Tetric EvoCeram may lead to sensitization. Tetric EvoCeram should not be used in such cases. To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (selectively apply a calcium hydroxide-based preparation in areas close to the pulp).

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Interactions

Substances containing eugenol/clove oil may inhibit the polymerization of composite materials. Consequently, the application of such materials together with Tetric EvoCeram must be avoided. Cationic mouth rinses and plaque disclosing agents as well as chlorhexidine may cause discolouration.

Application

1. Shade determination

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist.

2. Isolation

Adequate relative or absolute isolation using suitable auxiliaries, such as OptraGate® or OptraDam® Plus, is required.

3. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e., by preserving as much of the tooth structure as possible. Do not prepare sharp internal edges and angles or additional undercuts in caries-free areas. The dimensions of the cavity are generally determined by the extent of the caries or the size of the old restoration. Bevel enamel margins of anterior teeth. In the posterior region, only the sharp enamel edges should be rounded (finishing diamonds, 25–40 µm). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes using rubber cups or rotary brushes. Subsequently, remove all residue in the cavity with water spray and dry with water- and oil-free air.

4. Pulp protection / Base

Do not apply a base material when using an enamel/dentin bonding agent. Selectively coat only very deep areas close to the pulp with a calcium hydroxide material (e.g. ApexCal®) and subsequently cover with a pressure-resistant cement (e.g. a glass-ionomer cement, such as Vivaglass® Liner). Do not cover the other cavity walls so that they can be used to support the bond with an enamel/dentin adhesive.

5. Apply matrix / interdental wedge

Use a wrap around matrix for cavities affecting the proximal area or a sectional matrix and wedge it.

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the Instructions for Use of the product in use. Ivoclar Vivadent recommends using the universal bonding agent Adhese® Universal, a single-component, light-cured adhesive for direct and indirect bonding procedures that features compatibility with all etching techniques.

7. Application of Tetric Evo Ceram

- For an optimum result, apply Tetric EvoCeram layers of max. 2 mm (or 1.5 mm Dentin shades, Bleach XL) and adapt with a suitable instrument (e.g. OptraSculpt).
- Sufficient exposure to the curing light prevents incomplete polymerization. For the recommendations regarding exposure time and light intensity see table 1.
- When using a metal matrix, additionally polymerize the composite material from the buccal or the lingual/palatal aspect after removing the matrix.
- A flowable composite (e.g. Tetric EvoFlow® or Tetric EvoFlow® Bulk Fill) can optionally be used. This optional layer must be cured separately (please refer to the respective Instructions for Use).

8. Finishing / Checking the occlusion / Polishing

Remove excess material with diamond or tungsten carbide finishers after polymerization. Check the occlusion and articulation and apply appropriate grinding corrections to prevent premature contacts or undesired articulation paths on the surface of the restorations. Use silicone polishers (e.g. OptraPol®) as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Additional notes

1. Tetric EvoCeram can be used in combination with Tetric EvoFlow Bulk Fill, Tetric EvoCeram Bulk Fill and Tetric EvoFlow.
2. In the case of repairs, additional Tetric EvoCeram can be directly applied to the polymerized material. If the Tetric EvoCeram restoration has already been polished, it must first be roughened and wetted with Adhese Universal before fresh Tetric EvoCeram is applied.
3. Tetric EvoCeram should have at least room temperature when applied. Cool temperatures render the material difficult to extrude.
4.  For single use only. If Tetric EvoCeram is applied from the Cavifil directly into the oral cavity of the patient, the Cavifil or the application cannula of the syringe should only be used for one patient due to hygienic reasons (prevention of cross-contamination between patients).
5. Syringes or Cavifils should not be disinfected with oxidizing disinfection agents.

Warning

Avoid contact of unpolymerized Tetric EvoCeram with skin, mucous membrane and eyes. Unpolymerized Tetric EvoCeram may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- Do not use Tetric EvoCeram after the indicated date of expiration.
- Storage temperature 2–28 °C (36–82 °F)
- Close syringes/Cavifils immediately after usage. Exposure to light leads to premature polymerization.
- Date of expiration: see information on Cavifils, syringes and packages.

Keep out of the reach of children!

For use in dentistry only!

The material has been developed solely for use in dentistry and should be used strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

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Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.