# Flow-It™ ALC

Flowable Composite

English

INSTRUCTIONS FOR USE

Français

French

MODE D'EMPLOI

Deutsch

German

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Italiano

Italian

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Español

Spanish

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Čeština

Czech

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English

INSTRUCTIONS FOR USE

## Indications:

- Flow-It ALC is a flowable visible light curing hybrid composite restorative material indicated for all Class I, II, III, IV, V restorations

## 1. Preparation

Prepare tooth as usual.

## 2. Etch and Apply Dentin Adhesive

Use Bond-1™ Primer Adhesive or other adhesive system as indicated by the instructions.

## 3. Place Flow-It ALC Flowable Composite

Apply Flow-It ALC Flowable Composite directly with the bent needle tip provided in the kit. Be sure to cover preparation site thoroughly with material. Do not build increments of more than 2mm thickness for proper light curing.\* Pull back on syringe plunger after use to prevent material seepage.

## 4. Light Curing

- With a high speed light:* Cure initial increments for 3 seconds, and the final increment for 6 seconds.
- With a standard halogen light:* Cure initial increments for 10 seconds, and the final increment for 20 seconds.

## 5. Finishing and Polishing

Use conventional method and materials for shaping. Use non-abrasive polishing disc with polishing paste to achieve the maximum shine.

**Optional:** For maximum surface sealing following polish, apply Protect-It™ Composite Sealant to seal the restoration surface and reduce wear.

\*Note: Due to their opacity, build up **Universal Opaque, White Opaque, Gingival Light and Gingival Dark** shades in increments of no more than 1mm thickness to ensure effective polymerization.

**NOTE:** Tips are for single patient use only, to prevent cross contamination between patients.

**CAUTION:** Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

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Store at ambient temperature.

## LIMITED WARRANTY —

Pentron's technical advice, whether verbal or in writing, is designed to assist dentists in using Pentron's product. The dentist assumes all risk and liability for damages arising out of the improper use of Pentron's product.

In the event of a defect in material or workmanship, Pentron will replace the defective product or part thereof, or reimburse the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Pentron.

## PENTRON

Manufactured by Pentron Clinical

1717 West Collins Ave. Orange, CA 92867 USA

800-551-0283 / 714-516-7557 / [www.pentron.com](http://www.pentron.com)

**EC REP** European Union Representative

SpofaDental a.s.

Markova 238

CZ-506 01 Jičín

Phone: +420 493 583 111