GC Fuji IX

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GC Fuji IX adalah bahan restoratif glass-ionomer (GIC) yang dirancang untuk teknik menambal gigi tanpa sakit/trauma (Atraumatic Restorative Treatment). The cavity preparation technique is quite different from normal operative dentistry and the application technique has been modified in order to allow the operator to achieve a high level of success with limited facilities.

ADVANTAGES

- Simple reliable technique for all operators.
- Chemical adhesion to dentin and enamel.
- Continuous fluoride release.
- Adequate placement time and rapid set.
- Excellent biocompatibility.
- Minimal water sensitivity
- High compressive strength
- Proven track record

RECOMMENDED INDICATIONS

- Class I and small Class II restorations in deciduous teeth.
- Class I and small Class II restorations in permanent teeth.
- Intermediate restorative and base material for heavy stress situation in Class I and Class II cavities using sandwich laminate technique.
- Class V and root surface restorations.
- Core build-up.
- Contra-indicating.

CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

PROPERTIES ON USAGE

Powder/Liquid ratio (g/g)	3.6/1.0
Mixing time (sec.)	25"-30"
Working time (min., sec)*	2'00"
Net setting time (min., sec)*	2'20"

Test conditions: Temperature (23±1-1°C), Relative Humidity (50±10%)

DIRECTIONS FOR USE



- Cavity Preparation
 - Using gloved fingers, isolate the tooth with cotton rolls, thus keeping the treatment area free from saliva..
 - Remove loose food debris and plaque with an explorer from the pits and fissures.
 - Wash the pits and fissures with a wet cotton pellet. The extent of the carious lesion will then become visible.
 - If the opening in the enamel is too small, widen the entrance by placing one edge of a Hatchet or Enamel Access Cutter in the entrance and turn the instrument around, making a scraping movement. Small particles of enamel will chip off. Remove the pieces of enamel by moving a wet cotton pellet around in the cavity. Dry the cavity with one or two pellets before continuing the procedure.
 - Carious dentine can now be removed. Be sure to wet the cavity now and then for easy removal of the loosened tooth material. It may be necessary to widen the entrance to the cavity further after the initial removal of carious dentin. Remove carious dentine with hand held excavators in a scooping movement, starting at the dentine-enamel junction and ending at the floor of the cavity (Fig. 1). The lower layers of carious dentin close to the pulp can be left behind.
 - Clean the cavity with a wet cotton pellet and dry thereafter with a dry pellet (Fig. 2).
 - Ensure that the adjacent pits and fissures are free from debris.
 - Make sure that all the demineralized enamel at the cavity opening is removed to ensure the tight sealing.
- Conditioning the cavity
 - Use the nozzle is clean and dry before dispensing the liquid.
 - Turn vertically and hold the tip about 5cm above the mixing pad, steady your hand and squeeze the bottle slowly and gently (Fig. 3).
 - Dispense 1 drop onto the corner of the mixing pad and use this for conditioning. Immediately dispense the second drop in the centre of the pad to be used later for mixing with the powder. Wipe the nozzle of the liquid bottle and replace the cap.
 - Take a clean cotton pellet, dip into the liquid reserved for conditioning and apply to the entire prepared cavity surface and adjacent pits and fissures for 10-15 seconds.
 - Make sure the pellet touches the cavity walls and use pellets appropriate to the size of the cavity.
 - Take a clean cotton pellet, dip in water, and wash the cavity 2 or 3 times for 5 seconds each.
 - Dry the cavity with a fresh cotton pellet. Surfaces should look shiny. The cavity and adjacent pits and fissures are now ready for the GC Fuji IX restorative material.
- Note:
 - Keep preparation dry and free from saliva and/or blood. Ensure proper isolation
- Dispensing
 - Powder to liquid ratio

The standard powder to liquid ratio is 3.6g/1.0g. This consistency can be obtained with one level spoonful of powder and one drop of liquid.
 - Dispensing powder

Tap the bottle on a hard surface to ensure there are no air pockets in the powder. Overfill the powder measuring spoon and level air with the powder with the collar at the top of the bottle and carry it into the mixing pad.
- Dispensing liquid
 - See directions in 2. above.
- Mixing

Using the plastic spatula provided, divide the powder into two equal parts. Spread the liquid across the mixing pad and mix the first half with the whole liquid for 10 seconds. Add the second part and mix for 15-seconds to obtain a homogeneous mixture (total time not more than 30 seconds) (Fig. 4).
- Placement
 - Insert the mixture into the cavity using the flat end of a carving instrument (Fig. 5).
 - Push the material with the round end of a medium sized excavator into the corner(s) using slight pressure making sure not to include any air bubbles.
 - Insert a second portion of the material and press into place with the round end of a large sized excavator.
 - Spread additional material on the occlusal surface and cover all the fissures adjacent to the cavity.
 - Make sure to complete this procedure while the material still has a glossy surface.
 - As soon as GC Fuji IX starts to lose the glossy surface rub some GC COCOA BUTTER or petroleum jelly over your index finger, place the finger over the tooth surface and apply finger pressure (Fig. 6). This will further compress the material into the cavity. The finger can be removed sideways from the tooth surface after about 30 seconds.
 - Remove visible excess material with the carver end of the applicator instrument
 - Check the occlusion
 - Wait until the material loses the shiny surface and adjust the occlusion with a medium sized excavator
 - Use the round end of an excavator to achieve a smooth material-tooth surface interface (Fig. 7).
 - Protect the restoration with a thin layer of GC COCOA BUTTER or petroleum jelly
 - Remove the cotton rolls used for isolation
 - Instruct the patient not to eat for at least 1 hour.

BRIEF REVIEW

			
Caries in occlusal surface of molar	Circular cavity preparations of the excavator	Cavity prepared using the ART technique	1.Cavity filled 2.Fissures filled and sealed

Note:

- Avoid use of GC Fuji IX with patients who have a history of allergy to glass ionomer cements.
- Do not use GC Fuji IX with patients who show an allergy to the material. In case of allergy, immediately stop using it and advise them to consult a physician.
- The operator who has a history of allergy to glass ionomer cements should not handle GC Fuji IX. In case of allergy, he should immediately stop using it and consult a physician.
- Do not allow the mixture to contact the oral tissues or skin. In case of contact, immediately remove the material with absorbent cotton soaked in alcohol and rinse with water.
- Avoid contact of the mixture with eyes. In case of contact, immediately flush with water and seek medical attention
- Clean the nozzle of the liquid bottle occasionally with a damp cotton pad to avoid clogging.
- If pulp capping is required, a calcium hydroxide such as GC Nu-Cap can be used to protect pulp.
- After use, tightly close both liquid and powder bottles to prevent exposure to moisture.
- Do not mix the powder or liquid of GC Fuji IX with any other glass ionomer product.
- The material is designed for use at temperature of 21-25°C (69-87.7°F). This product is for use only by a dental professional.
- This product should be used only for the applications described in the instructions.

SHADE

Universal shade

PACKAGE

GC Fuji IX 1-1 pack:
1 x 10g powder, 1 x 6g (4.8mL) liquid, 2 x 2g GC COCOA BUTTER, plastic spatula, powder scoop, and a mixing pad (#22)

GC Fuji IX 1-1 mini pack:
1 x 5g powder, 1 x 4g (3.2mL) liquid, 1 x 2g GC COCOA BUTTER, powder scoop

STORAGE

Store in a cool and dark place (4-25°C / 39.2-77.0°F). Shelf life
GC Fuji IX powder and liquid: 3 years from date of manufacture
GC COCOA BUTTER: 5 years from date of manufacture

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