ENGLISH

General Information

3M™ ESPE™ Ketac™ Nano Light-Curing Glass Ionomer Restorative, is a light curable, esthetic, and radiopaque (ISO standard 9917-2 equivalent 2.6 mm of aluminum), two part paste/paste, resin modified glass ionomer restorative, and primer system. Ketac Nano restorative provides the major benefits typical of other glass ionomer materials, along with the additional benefits of improved esthetics, wear resistance and polish. Ketac Nano restorative consists of a two part paste, packaged in a direct delivery Quick Mix Capsule. The two part paste offers greater convenience over the traditional powder/liquid hand-mix restorative systems.

Ketac Nano restorative is composed of unique fillers that provide benefits associated with conventional glass ionomers, resin modified glass ionomers and nanofilled composites. The filler content of this system is a combination of two types of surface treated silica/zirconia nanofillers approximately 5-25 nm and nanoclusters approximately 1.0 to 1.6 microns along with a radiopaque fluoroaluminosilicate glass of an average particle size of approximately 1 micron. The inorganic filler loading is about 69% by weight or (56% by volume).

Ketac Nano restorative **must** be used with Ketac Nano Glass Ionomer Primer, a one part, visible light cure cavity primer. Its function is to adequately wet the bonding surfaces to facilitate adhesion of the Ketac Nano resin modified glass ionomer restorative. In use, the primer is dispensed, applied, air dried and light cured. Adequately air drying and then light curing the primer separately before placement of the Ketac Nano resin modified glass ionomer is required for adhesion of the glass ionomer to tooth structure.

Indications

Ketac Nano restorative system is indicated for:

- Primary teeth restorations
- Small Class I restorations
- Class III and V restorations
- Transitional restorations
- Filling defects and undercuts
- Laminate/Sandwich technique
- Small core build-ups where at least 50% of coronal tooth structure is remaining for support

Generally coronal tooth structure that is exposed to high stress factors, i.e. cusps, should not be restored with any glass ionomer type restorative materials.

Contraindications

The Ketac Nano system **is not recommended for direct pulp capping**. Cover areas in close proximity of a pulpal exposure with a small amount of hard setting calcium hydroxide material, or a resin modified glass ionomer liner (e.g. 3M™ ESPE™ Vitrebond™ Plus Light Cure Glass Ionomer Liner/Base).

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Precautionary Information for Patients:

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

Precautionary Information for Dental Personnel:

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed.

Warning: The Ketac Nano Quick Mix Capsule must be used with a 3M™ ESPE™ Aplicap™ applier, or an applier of identical design. Use of an applier of a different design (e.g. those designed for jumbo tubes) can cause serious injury because capsule may eject forcefully from applier.

3M ESPE MSDS can be obtained from www.3MESPE.com or contact your local subsidiary.

Instructions for Use

1. Shade selection. For esthetic restorations, select the desired shade using the Ketac Nano Light Curing Glass Ionomer Restorative shade guide. Refer to label on each individual foil package for shade identification.

2. Cavity Preparation. Remove caries. Rinse and slightly dry cavity, do not desiccate. Ketac Nano restorative is not recommended for direct pulp capping. Apply a liner if pulp tissue is near exposure conditions.

3. Primer. Ketac Nano Glass Ionomer Primer **must** be applied to the preparation before restoring tooth with Ketac Nano restorative. Core buildups with missing cusps may require placement of pins for retention. For core buildups with pins, apply primer to pins as well.

3.1 Dispense the Ketac Nano primer into a well.

3.2 Using a fiber tip, apply primer for 15 seconds to prepared semi-dry enamel and dentin surfaces. Replenish primer as needed to assure surfaces are kept wet with primer for the recommended application time.

3.3 Dry the primer using an air syringe for 10 seconds. **Do not rinse.** After drying, the primed surfaces will remain shiny in appearance.

3.4 Light cure the primed surfaces for 10 seconds using a 3M ESPE curing unit or other dental visible light curing unit of comparable intensity.

3.5 The light cured surfaces will appear shiny.

Note: By adequately drying and separately light curing the primer, maximum adhesion of the glass ionomer to tooth structure is obtained. The primer is light sensitive and contains photoinitiators and water. Minimize ambient light exposure and evaporation by dispensing just prior to use and replacing vial cap immediately after dispensing.

4. Dispensing Quick Mix Capsule

The Quick Mix Capsule will deliver approximately 0.3 grams of restorative material.

Note: A fail safe feature has been designed into the Quick Mix Capsule to prevent the operator from over pressurizing it and causing possible injury. The following actions will result in the leakage of paste at the base of the orange nozzle.

- Failure to lift the nozzle into the full open position
- Dispensing paste too fast from the capsule
- Attempting to dispense paste beyond the 90 second extrusion time.

4.1 Just prior to use, remove Quick Mix Capsule from foil package.

4.2 When ready to dispense Ketac Nano restorative into preparation, lift orange mix tip until it is in a straight line with the capsule. Do not force beyond stop. Once the nozzle is swung open and activated **do not** reclose as this may cause a capsule failure.

4.3 Place capsule into the applier gun such that the capsule holder engages the groove at the plunger end of the capsule. Press the capsule down into the holder as far as it will go. **Do not** engage the applier plunger into capsule until ready to use.

4.4 To dispense Ketac Nano squeeze handle **slowly** to extrude a small amount of material approximately 2-3 mm in diameter outside the mouth to verify capsule function. Discard this material. Time to dispense paste from Quick Mix Capsule is 90 seconds. Exceeding the 90 second time may effect properties of Ketac Nano restorative or cause capsule failure.

4.5 Quick Mix Capsule is now ready to dispense material directly into the preparation. The capsule can be rotated inside the holder to achieve an optimal position for application. Keep tip immersed in material to minimize air entrapment.

4.6 **Attention:** If placing material incrementally, additional material extruded from the Quick Mix Capsule must be delivered into the preparation within the 90 second dispensing time.

5. Placement.

Working time of Ketac Nano restorative is 2.5 minutes beginning with the extrusion of material into the preparation.

Quick Mix Capsule	2.5 minutes working time
-------------------	--------------------------

In a semi-dry to dry field incrementally syringe or place material in a depth of 2 mm or less.

Handling tips

- Dental instruments used for shaping and contouring can be moistened with Ketac Nano Primer to prevent glass ionomer restorative from adhering to them. Do not saturate the restoration with primer.

- Allowing Ketac Nano restorative to rest for approximately 60 seconds after placement in preparation permits the paste to become firmer and less tacky. After this waiting period the material may then be shaped and contoured with limited or no adhesion to the placement instrument.

- Note:** Waiting time should not exceed the total working time.

5.1 **Cavities.** Syringe the mixed Ketac Nano restorative into the cavity keeping the syringe tip immersed in the material to minimize air entrapment. Material will become firmer, shortly after placement. Shape and contour the restoration using an appropriate placement instrument.

5.2 **Small Core buildups.** Syringe the Ketac Nano restorative into undercut areas, around pins, around posts, and then fill the preparation in increments of 2 mm or less.

Note: Once cured the prepared Ketac Nano core buildup is compatible with conventional impressioning materials. The core buildup should be kept wet with saliva or lubricated to prevent bonding to chemical cure provisional. The glass ionomer core buildup will not bond with temporary cements.

5.3 Laminates/Sandwich Technique. Use when cavity design allows for a minimum composite restorative thickness of 2 mm on occlusal surface. Syringe material into the preparation if margins are located partially in dentin or aprismatic enamel as, for example, in deep Class II cavities. Syringe material in prepared cavity extending the glass ionomer restorative base no further than apical to the proximal contact point.

Note: After light curing refine by removing excess Ketac Nano primer and Ketac Nano restorative base material from enamel margins and cavity walls to be bonded subsequently with the adhesive/composite system.

6. Curing. Ketac Nano restorative will cure only by exposure to visible light. **The maximum depth of material for light curing should not exceed 2 mm.** Light cure the Ketac Nano restorative by exposing its entire surface area to a 3M ESPE visible light curing unit or other dental light curing unit of comparable intensity according to the chart.

3M ESPE Light	mW/cm²	Seconds	Shades
Elipar™ S10 LED	≥1200	20	All
Elipar™ FreeLight 2 LED	≥1000	20	All
Elipar™ 2500 Halogen	≥450	20	A1, A2, A3, B2
Elipar™ 2500 Halogen	≥450	30	A3.5, A4

7. Finishing. Immediately after curing, the Ketac Nano restoration can be contoured and polished using conventional finishing and polishing instruments, (e.g. Sof-Lex™ Finishing and Polishing System), under moist conditions. Do not allow surface to become desiccated.

Storage and Use

- This product is designed to be stored and used at room temperature. Ambient temperature routinely higher than 27°C/80°F, or lower than 2°C/35°F may reduce shelf life. See outer package for expiration date.

- Quick Mix Capsules removed from sealed foil packs must be used within 24 hours.

- Quick Mix Capsules are intended for single-use only and must be discarded after use.

- Primer and paste are light sensitive materials. Protect them from ambient light exposure by dispensing just prior to use and replacing caps immediately after dispensing.

- Disinfect primer bottle using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association, *Guidelines for Infection Control in Dental Health-Care Settings –2003* (Vol. 52; No. RR-17), Center for Disease Control and Prevention.

Disposal

See the Material Safety Data Sheet (available at www.3MESPE.com or through your local subsidiary for disposal information).

Customer information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Caution: U.S. Federal Law restricts the device to sale or use on the order of a dental professional.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the products for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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